

Introduction

Phosphoprotein Gel Stain is a high sensitive fluorescent stain designed for selectively detecting phosphoproteins in polyacrylamide gels. This stain contains a phos-tag™ group, which allows direct, in-gel detection of phosphate groups attached to tyrosine, serine, or threonine residues, without the need for antibodies or radioisotopes. The stain can be used with standard SDS-polyacrylamide gels or with 2-D gels.

Package Information

Components	C0118
Phosphoprotein Gel Stain	500 ml

Ex/Em: 550/580 nm

Storage

Store at 2-25°C, Protect from light

Sample Preparation

A delipidated and desalted sample is essential for adequate separation of proteins by electrophoresis and subsequent staining by Phosphoprotein Gel Stain.

- For a 150 µl sample (~150-300 µg of protein), add 600 µl of methanol and mix well by vortexing.
- Add 150 µl of chloroform and mix well by vortexing.
- Add 450 µl of ultrapure water and mix well by vortexing.
- Centrifuge at ~12,000 rpm for 5 min.
- Discard the upper phase, keeping the white precipitation disc that forms between the upper and lower phases.
- Add 450 µl of methanol and mix well by vortexing.
- Centrifuge at ~12,000 rpm for 5 minutes.
- Discard the supernatant and dry the pellet in a vacuum centrifuge for 10 minutes.
- Resuspend the pellet in standard 1× sample buffer for electrophoresis.

Staining Protocol

Note: The protocol is optimized for standard 1 mm thick, 8 cm × 8 cm SDS-PAGE minigels. Larger or thicker gels require additional volumes of reagents or longer incubation times.

- Run gel as usual according to your standard protocol.
- Fix gel with 100 ml of fix solution (50% methanol, 10% acetic acid), and agitate on an orbital shaker for 30 min. Repeat one more time with 100 ml fresh fix solution.
- Wash the gel in 100 ml of ultrapure water with gentle agitation for 10 minutes. Repeat this step twice, for a total of three washes.

Phosphoprotein Gel Stain

Cat. #: C0118 Size: 500 ml

- Stain the gel with enough Phosphoprotein Gel Stain (40~60 ml) to cover the gel, and agitate on an orbital shaker for 60-90 min.
- Destain the gel with Phosphoprotein Destain Solution (Cat. #: C00119) with gentle agitation for 30 minutes. Repeat this procedure two more times.
- Wash the gel twice with ultrapure water for 5 minutes per wash. If the background is high or irregular, the gel may be left in the second wash for 20-30 minutes and re-imaged.
- Image gel using recommended instruments and filter sets (see Table 1 for recommendations). A 300 nm UV transilluminator or a blue-light transilluminator can be also used for imaging. However, the sensitivity will be 10-fold lower.

Protocol Quick Reference

	Reagent	Protocol
Fix	50% methanol, 10% acetic acid	100 mL, 30 min
		100 mL, 30 min
Wash	Ultrapure water	100 mL, 10 min
		100 mL, 10 min
		100 mL, 10 min
Stain	Phosphoprotein Gel Stain	40-60 mL 60-90 minutes.
Destain	Phosphoprotein Destain Solution	60 mL, 30 min
		60 mL, 30 min
		60 mL, 30 min
Wash	Ultrapure water	100 mL, 5 min
		100 mL, 5 min

Staining the Gel for Total Protein

After staining with Phosphoprotein Gel Stain, the gel can be stained with a total-protein stain.

- Image the gel following staining with the first gel stain.
- Rinse the gel with ultrapure water for 5 minutes. Repeat one more time.
- Incubate gel with SuperLumin™ Protein Gel Stain (Cat. #: C00117) solution (40~60 ml). Microwave 45 seconds, and agitate on an orbital shaker for 15 min. Repeat microwave 45 seconds, and agitate on an orbital shaker for another 15 min.
- Wash gel with 100 ml wash solution (10% methanol, 7% acetic acid) for 30 min.
- Image gel with a 300 nm UV transilluminator, blue-light transilluminator or a laser scanner.

Table 1. Filters recommended for use with Phosphoprotein Gel Stain

Instrument	Manufacturer	Excitation Source	Emission Filter
Typhoon Trio+, Trio, 9200, 9210, 9400, 9410	Amersham Biosciences	532 nm laser	560 nm longpass
FluorImager	Amersham Biosciences	514 nm laser	570 nm bandpass
Molecular Imager FX	Bio-Rad Laboratories	532 nm laser	555 nm longpass
FLA-3000G, FLA-5100	Fuji Photo Film Co, Ltd	532 nm laser	580 nm longpass
ProXPRESS	PerkinElmer LifeSciences	540/25 nm	590/30 nm